AEROSOL



Micro Technologies, Big Ideas.

PRODUCT INFORMATION

Model 135—Mini-MOUDI™ Impactor

A 2 L/min precision cascade impactor for collecting size-fractionated particle samples for gravimetric and/or chemical analyses



INTRODUCTION

Model 135 cascade impactors are designed for sampling and collecting size-fractionated aerosol particle samples for gravimetric and/or chemical analyses. The various versions of the Model 135 impactors differ only by the number of stages, the last stage's cutpoint size, and the type of inlet.

DESCRIPTION

The Model 135 Mini-MOUDI impactors can have 6, 8 or 10 stages, operating at a sampling flow rate of 2 L/min for a variety of sampling applications. The Mini-MOUDI impactor has stages available with cut-point diameters of 0.056, 0.10, 0.18, 0.32, 0.56, 1.0, 1.8, 3.2, 5.6 and 10 μ m.

For a Model 135 stage, half of the area is the impaction plate for the stage above while the other half contains the nozzles for that stage. Subsequent stages alternate in orientation so that there is always an impaction surface below each set of nozzles. This arrangement allows for a very compact cascade impactor design (Figure 2).


The Model 135 impactors can accept a variety of sampling substrates, including aluminum foils, thin plastic films and membrane filters. The substrates are 37-mm in diameter with a semicircular hole in one half of the substrate (Figure 3).

The stages are designed to prevent cross-flow interference between adjacent nozzles. The result is sharp cut-size characteristics not available with other cascade impactors that are less well designed aerodynamically (Figure 3).

The Models 135-6 and 135-8 can be used with personal sampling pumps since they have a low pressure drop. The Model 135-10 requires a vacuum pump due to the larger pressure drop associated with the last two stages.



Figure 2. Impaction substrate



FEATURES

- 2-L/min sampling flow rate
- Up to 10 impactor stages with nominal cutpoints at 0.056, 0.10, 0.18, 0.32, 0.56, 1.0, 1.8, 3.2, 5.6 and 10 µm, plus a final filter
- 37-mm diameter semicircular impaction substrates
- 37-mm final filter
- Impactor is made of anodized aluminum for light weight, durability and nozzle dimensional stability
- Low pressure drop
- Sharp cut-off characteristics
- · Open Inlet or Inlet Cone
- Low inter-stage losses

APPLICATIONS

- Personal exposure particle monitoring
- Work place aerosol analysis
- · Cigarette and E-Cig smoke characterization
- Indoor residential aerosol studies

TO ORDER

Mini-MOUDI Impactors:

| Specify | Description | |
|---------|---|--|
| 135-6A | Mini-MOUDI Impactor, 6 Stage | |
| 135-6B | Mini-MOUDI Impactor, 6 Stage, Inlet Cone | |
| 135-8A | Mini-MOUDI Impactor, 8 Stage | |
| 135-8B | Mini-MOUDI Impactor, 8 Stage, Inlet Cone | |
| 135-10A | Mini-MOUDI Impactor, 10 Stage | |
| 135-10B | Mini-MOUDI Impactor, 10 Stage, Inlet Cone | |

Accessories:

| Specify | Description | | |
|--------------|--|--|--|
| 0135-75-5007 | Vacuum Pump, 135-6,135-8, 110V charger | | |
| 0135-75-5008 | Vacuum Pump, 135-6,135-8, 220V charger | | |
| 0135-01-0100 | Vacuum Pump, 135-10 , 110V | | |
| 0135-01-0101 | Vacuum Pump, 135-10 , 220V, EU | | |
| 0135-01-0102 | Vacuum Pump, 135-10 , 220V, UK | | |
| 0135-01-0014 | Al Foil Substrates, Pkg 500 | | |
| 0135-01-5203 | Glass Fiber Filters, 37-mm, Pkg 100 | | |

SPECIFICATIONS*

*Specifications are subject to change without notice

| | Model 135-6 | Model 135-8 | Model 135-10 |
|--|------------------------------|---|---|
| Impactor Stages | 6 | 8 | 10 |
| Pressure Drop, kPa (w/o final filter) | 2.5 | 10 | 40 |
| Cut-point diameters, µm | 0.56, 1.0, 1.8, 3.2, 5.6, 10 | 0.18, 0.32, 0.56, 1.0, 1.8, 3.2, 5.6,10 | 0.056, 0.10, 0.18, 0.32, 0.56, 1.0, 1.8, 3.2, 5.6, 10 |
| Dimensions (D x H) | 66 x 62 mm | 66 x 77 mm | 66 x 85 mm |
| Weight (total) | 0.225 kg (0.50 lb.) | 0.325 kg (0.72 lb.) | 0.415 kg (0.91 lb.) |
| Vacuum Pump | 4.8V, 2.15 A-h NiMH battery | 4.8V, 2.15 A-h NiMH battery | 115 or 230 VAC, 50-60Hz, 0.07 kW |

Environmental Operating Conditions

10-50°C (50-122°F)

0%-90% RH (non-condensing)

MSP Corporation

5910 Rice Creek Parkway, Suite 300 Shoreview, Minnesota 55126, U. S. A. Phone: 651.287.8100; Fax: 651.287.8140 sales@mspcorp.com; www.mspcorp.com

Copyright© 2017 MSP Corporation (MSP-PI-135 Rev. A – A4). All rights reserved.