

# DIM KITS FOR BSC APPLICATIONS

## FOR MODELS EBT731 BALOMETER<sup>®</sup>, PH731 PROHOOD<sup>™</sup>, AND 8380 ACCUBALANCE<sup>®</sup> CAPTURE HOODS

**Comply with NSF 49  
the Easy Way**

APPLICATION NOTE TSI-152

To certify Class II Bio-Safety Cabinets, just use the Direct In-Flow Measurement Hood. Designed by TSI for certifying professionals, this easy-to-use tool allows you to take accurate measurements fast and document them for NSF 49 compliance.

### Throw away the cardboard!

Expansion panels slide out to cover openings from 2 to 6 feet wide. Choose the BSC hood and frame kit that best suits your work.

BSC Hood and Frame Kit Models:

- 801204: 8-in x 21-in (200 mm x 530 mm)
- 801205: 10-in x 21-in (255 mm x 530 mm)



BSC Hood kits come with an adjustable tripod and mounting instructions for easy setup.

### Simply attach the BSC hood kit to the instrument base and take full advantage of all its features:

- Ergonomically designed base and meter
- Extremely light weight
- Large multi-line back-lit display
- Wheeled luggage-style carrying case

### Save valuable time:

- Make and log readings in 2 to 8 seconds
- Data log up to 26,500 samples
- Download to PC with included software

### Save money:

- Low purchase price
- Low service costs
- Modular design allows you to buy the basic kit and add components as needed and budget allows





UNDERSTANDING, ACCELERATED

**TSI Incorporated** – Visit our website [www.tsi.com](http://www.tsi.com) for more information.

<b>USA</b>	<b>Tel:</b> +1 800 874 2811	<b>India</b>	<b>Tel:</b> +91 80 67877200
<b>UK</b>	<b>Tel:</b> +44 149 4 459200	<b>China</b>	<b>Tel:</b> +86 10 8251 6588
<b>France</b>	<b>Tel:</b> +33 4 91 11 87 64	<b>Singapore</b>	<b>Tel:</b> +65 6595 6388
<b>Germany</b>	<b>Tel:</b> +49 241 523030		