

Rotating Disk Thermodiluter

Model 379020

Thermal Conditioner

Model 379030

Standard replacement parts for Model 379020 Rotating Disk Thermodiluter consist of:

| <u>Part Number</u> | <u>Part Description</u> |
|---------------------------|--|
| 379025 | Housing for 379020 Diluter Head |
| 1137001 | Particle Filter for 379020/21/30/31 |
| 1137003 | Plastic tube for sample gas (3 m) |
| 1137004 | Silicon tube for undiluted gas (3 m) |
| 1137005 | Peristaltic Pump 60L/h 24 VDC |
| 1137010 | Roller for peristaltic pump |
| 1137006 | Tubing for new peristaltic pump (98019) |
| 1137002 | Tubes for old peristaltic pump |
| 1137011 | O-rings for 3m connection tubes (qty. 6) |
| 1137009 | 8 cavity disk w/calibration |
| 1137008 | 10 cavity disk w/calibration |
| 1137007 | Dilution block w/calibration |
| 1132012 | 8 or 10 cavity disk calibration |

Standard replacement parts for Model 379030 Thermal Conditioner consist of:

| <u>Part Number</u> | <u>Part Description</u> |
|---------------------------|-------------------------------------|
| 1137001 | Particle Filter for 379020/21/30/31 |
| 1137301 | Rotary Valve Pump 21 L/min, 24 VDC |
| 1137302 | Rotary Valve Pump 3.1 L/min, 12 VDC |
| 1137303 | 379030 filter 0.02um/99.9999% |

For pricing information, or to place an order you can 1) send the necessary information to answers@tsi.com or 2) contact your local TSI representative (see contact information below) and speak to one of our customer specialists. Please include your contact information (telephone number and mailing address) with all requests.

TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

| | | | |
|----------------|-----------------------|--|---|
| USA | Tel: +1 800 874 2811 | E-mail: particle@tsi.com | Website: www.tsi.com |
| UK | Tel: +44 149 4 459200 | E-mail: tsiuk@tsi.com | Website: www.tsiinc.co.uk |
| France | Tel: +33 491 95 21 90 | E-mail: tsifrance@tsi.com | Website: www.tsiinc.fr |
| Germany | Tel: +49 241 523030 | E-mail: tsigmbh@tsi.com | Website: www.tsiinc.de |
| India | Tel: +91 80 41132470 | E-mail: tsi-india@tsi.com | |
| China | Tel: +86 10 8260 1595 | E-mail: tsibeijing@tsi.com | |

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.



TRUST. SCIENCE. INNOVATION.