

QUICK START GUIDE

STANDARD ACCESSORY KIT

5200 SERIES



P/N 6011692, REVISION A
APRIL 2019

TSI and TSI logo are registered trademarks of TSI Incorporated.
FLO-Sight is a trademark of TSI Incorporated.



For further information/customer support, please contact:

U.S. & INTERNATIONAL

Sales and Customer Service:

(800) 680-1220 / +1(651) 490-2860

Fax:

+1(651) 490-3824

Email:

answers@tsi.com

TSI Instruments Ltd. (UK)

Sales and Customer Service:

+44 (0) 1494 459200

Fax:

+44 (0) 1494 459700



5200 Series Standard Accessory Kit (52000)

Kit Components






TSI Part Number	Description	Picture
5200-025IN	5200 Series 1/4-inch Tube Ends (2)	
5200-10MM	5200 Series 10-mm Tube Ends (2)	
5200-IF-0375IN	5200 Series Inlet Filter 3/8-inch—Low Pressure	
5200-IF-025IN	5200 Series Inlet Filter 1/4-inch—Low Pressure	
5000-USBC-A	5000 Series USB-C to USB-A Power and Communications Cable (6 ft.)	

TSI Part Number	Description	Picture
5000-UPS	5000 Series Universal Power Supply	
5000-HUB	5000 Series USB Hub Cable	

5200 Series Optional Accessories

TSI Part Number	Description	Picture
5200-EC	5200 Series End Caps (2)	
5200-C	5200 Series Collars (2)	

5000 Series Optional Accessories


TSI Part Number	Description	Picture
5000-TTS	5000 Series Tabletop Stand	 <p>A tan-colored metal tabletop stand with a flat base and a vertical support. Three screws are shown below the stand.</p>
5000-RS232	5000 Series USB-RS232 Cable with Null Modem Converter	 <p>A white cable with a blue USB-A connector on one end and a white RS-232 connector on the other. A separate white null modem converter is shown next to it.</p>
5000-USB-LOCK	5000 Series USB-C to USB-A Cable with Screw Lock (Includes reversible head screwdriver)	 <p>A black cable with a USB-C connector on one end and a USB-A connector on the other. A small white reversible head screwdriver is shown next to it.</p>
5000-WMK	5000 Series Wall Mounting Kit	 <p>A metal wall mounting kit consisting of a main plate, a smaller plate, and several screws.</p>
5000-WMDRK	5000 Series Wall Mounting DIN Rail Kit	 <p>A metal wall mounting kit for DIN rails, including a main plate, two DIN rail brackets, and several screws.</p>
5000-DCC	5000 Series Deluxe Carrying Case	 <p>A blue and grey soft carrying case with a handle and a TSI logo on the front.</p>
5000-PC-ADV	5000 Series PC Software (Advanced version)	See following page for more details

FLO-Sight PC Software

Get the most out of your 5000 Series Gas Flow Multi-Meter with TSI's FLO-Sight™ PC Software.

Couple the real-time data visualization of FLO-Sight PC software with the fast response rate of our 5000 Series Flow Meters to gain valuable insights and truly understand your flows.

See your flows as you have never seen them before with FLO-Sight software.



FLO-Sight software has these basic features:

- Real-time display of all measurement parameters
- Min, max, and average readings
- Advanced meter configuration
- Volume measurements with user-set begin/end triggers (5210, 5220, 5230 models only)

FLO-Sight software advanced features include:

- Real-time graphing
- Data logging, download, and file management
- Alerts set for any parameter
- Multiple meter operation

A free trial of FLO-Sight software's advanced features is included with the Series 5200 Multi-Meter. The advanced version can be purchased through the software itself or at www.tsi.com.

For more information, contact: answers@tsi.com.



UNDERSTANDING, ACCELERATED

TSI Incorporated – Visit our website www.tsi.com for more information.

USA Tel: +1 800 680 1220

India Tel: +91 80 67877200

UK Tel: +44 149 4 459200

China Tel: +86 10 8219 7688

France Tel: +33 1 41 19 21 99

Singapore Tel: +65 6595 6388

Germany Tel: +49 241 523030

