

# AeroTrak<sup>®</sup> Remote Particle Counters With Pump

Models 6310, 6510 & 6510-VHP

Model 6310

AFROTRAK

4

VS.

1111

AeroTrak® Models 6310, 6510 and 6510-VHP are independent monitoring systems condensed into an easy-to-clean 316L stainless steel enclosure.

The integrated clean pump eliminates the need for an external vacuum source. Four standard analog inputs help to reduce installation and integration costs. Vaporized Hydrogen Peroxide (VHP) resistance and sample inlet capping auto-detect means this instrument is ideal for use in any GMP controlled area where continuous particulate and environmental monitoring is required.

Data is output via integrated Ethernet (TCP/IP), serial Modbus RTU communications and 4-20 mA analog outputs. Removing the unique hot-swap combined particle sensor and blower assembly from the mounted enclosure requires no tools, making it quick and easy to maintain. For ease of use, the internal pump automatically stops when the sample inlet probe is capped. TSI's optional FMS Buddy Automatic Hot-Standby Monitoring Software adds increased levels of assurance by maintaining full system functionality in the event of a monitoring computer failure.

The 3,000 sample record data buffer provides data redundancy, assuring full compliance and peace of mind during an unexpected network failure. Backed by TSI's comprehensive warranty and commitment for superior service and support, these remote particle counters lead the way with such a unique combination of features and benefits.

### Monitoring Applications:

- GMP Controlled Areas
- Isolators
  - Restricted Access Barrier Systems (RABS)
  - Cleanrooms
  - Hospital Pharmacies
  - Laboratories
  - Industrial
  - Aerospace

#### **Specifications**

# AeroTrak® Remote Particle Counters With Pump

Models 6310, 6510 & 6510-VHP

#### Easy to Install and Maintain

- · Eliminates the need for centralized vacuum system
- Internal blower has no carbon vanes and low heat gain Standard four 4-20 mA analog inputs reduce environmental
- sensor integration and installation costs Hot-swap optical sensor and pump slide out as a single
- assembly for fast and easy maintenance, no tools required Sample inlet capping auto-detect automatically stops the blower
- and prevents accidental damage
- Model 6510-VHP is resistant to VHP exposure, enabling use where VHP surface disinfection is utilized
- Easy-to-clean passivated 316L stainless steel enclosure, sample probe and cap, which are resistant to cleaning solutions
- Ethernet (TCP/IP) and serial Modbus RTU data output enables seamless integration into most software, including TSI's FMS Software

#### **Meets Current GMP Guidelines** and International Standards

- 1 CFM (28.3 L/min) flow rate, up to four channels of particle data including 0.5  $\mu m$  and 5  $\mu m$  per EU and US FDA Aseptic **Processing Guidelines**
- ISO 21501-4 compliant calibration delivers consistent instrument counting and sizing performance
- FMS 5 Monitoring Software integration enables 21 CFR Part 11-compliant data collection



Hot Swap Optical Sensor and Pump

Specifications are subject to change without notice.

AeroTrak, TSI and the TSI logo are registered trademarks in the United States and may be protected under other country's trademark registrations.

## 6310

```
Size Range
Channel Sizes
Counting Efficiency
```

0.3 to 25 µm 0.3, 0.5, 0.7, 1.0 µm 0.3, 0.5, 1.0, 5.0 µm 50% at 0.3  $\mu$ m; 100% for particles > 0.45  $\mu$ m (per ISO 21501-4 and JIS)

#### 6510 & 6510-VHP

Size Range **Channel Sizes Counting Efficiency**  0.5 to 25 µm 0.5, 0.7, 1.0, 5.0 µm 50% at 0.5 µm; 100% for particles > 0.75 µm (per ISO 21501-4 and JIS)

#### 6310, 6510 & 6510-VHP

821,000 particles / ft <sup>3</sup> (29,000,000/m <sup>3</sup> ) @ 10% coincidence loss
< 15% @ 0.5 µm (per ISO 21501-4)
<1 count / 5 minutes (<2 particles /ft <sup>3</sup> )
(per ISO 21501-4 and JIS)
1 CFM (28.3 LPM) with $\pm$ 5% accuracy
Blower*, -non particle generator, auto stops if
flow blockage detected due inlet capping – HEPA filtered exhaust
Long Life Laser diode
Ethernet (TCP/IP) and serial Modbus
output 4-20 mA Analog Outputs: 2 channel
log/linear particle counts, one channel
unit status
3,000 sample record buffer in the event of
network failure
4-20 mA, supports up to 4 environmental
sensors (loop powered)
Normal open dry contact rated 0 to
60 V AC/DC 1 Amp
Four LEDs: power, flow, status and sampling
FMS 5 Facility Monitoring Software, enables
21 CFR Part 11 Compliance
Via Ethernet (TCP/IP) Web interface
10.5 in. x 10 in. x 6.1 in.
(26.7 cm x 25.4 cm x 15.5 cm)
14.3 lb (6.5 kg)
24 VDC, 3 A
316L Passivated Stainless Steel
CE, JIS B9921, ISO 21501-4
Two Years, 1 Year (Model 6510-VHP)
36° to 104°F (2° to 40°C) / 20% to 95%
non-condensing
14° to 122°F (-10° to 50°C) / 20% to 98%
non-condensing
316L Stainless Steel Isokinetic Sample
Probe Inlet, Probe Cap and Chain, 316L
Stainless Steel Inlet, Flexible Sample Tubing,
Operating Manual on CD
24 VDC Power Supply, Zero Filter



TSI Incorporated - Visit our website www.tsi.com for more information.

USA	Tel: +1 800 874 2811	India	Tel: +91 80 67877200
UK	Tel: +44 149 4 459200	China	Tel: +86 10 8219 7688
France	Tel: +33 1 41 19 21 99	Singapore	Tel: +65 6595 6388
Germany	Tel: +49 241 523030		