



CERTIFICATE OF APPROVAL

This is to certify that the Quality Management System of:

TSI Instruments Singapore Pte Ltd
150 Kampong Ampat
#05-05 KA Centre
Singapore 368324
Republic of Singapore

has been approved by Lloyd's Register Quality Assurance (Singapore)
to the following Quality Management System Standards:

ISO 9001:2008

The Quality Management System is applicable to:

**Sale and service of measurement instruments and
accessories for the cleanroom industry.**

Approval
Certificate No: SNG 6007090

Original Approval: 24 August 2009

Current Certificate: 24 August 2009

Certificate Expiry: 23 August 2012



Issued by: Lloyd's Register Quality Assurance Limited

This document is subject to the provision on the reverse
71 Fenchurch Street, London EC3M 4BS United Kingdom. Registration number 1879370
This approval is carried out in accordance with the LRQA assessment and certification procedures and monitored by LRQA.
The use of the UKAS Accreditation Mark indicates Accreditation in respect of those activities covered by the Accreditation Certificate Number 001
Macro Revision 13